

# APPLICATION SHEET


Polymers – Pharmaceuticals – Photo DSC 204 **F1 Phoenix**<sup>®</sup>

## Methacrylates Used as Dental Restorative Material

### Introduction

In dental applications, light-curing dental composites are used as restoratives (fillings) or veneering materials. The materials are generally based on methacrylate systems such as bis-glycoldimethacrylate (bis-GMA) or urethandi-methacrylate (UDMA). Additional monomers are used as

diluents or to guarantee the cross-linking abilities of the resin. Additionally, inorganic fillers up to 80 weight% are used to improve the mechanical properties and reduce polymerization shrinkage. The examples for dental restorative composites measured here are Charisma (Heraeus Kulzer), Adamant and EvoCeram (Icovlar-Vicadent) and Admiro (Voco).



### Test Conditions

Temperature range: 51°C isothermal  
Heating/cooling rates: 0 K/min  
Atmosphere: Nitrogen (50 ml/min)  
Sample mass: approx. 25.5 mg  
Crucible: open Al  
UV device: Delolux 04  
Radiation time: 60 s

### Test Results

All materials showed a strong exothermic reaction after the UV light is activated. After the initial curing process, the UV system was switched on a second time. From the difference signal between the first and second illumination process, the reaction enthalpy was calculated as area below the peak. As a result of the different composition and filler content, different reaction enthalpies were found for the different restorative materials (EvoCeram 141.6 J/g, Charisma 128.4 J/g and Adamant 105.5 J/g).